


REGEN+U



SHAPER DATASHEET
2020

PRODUCT NAME	SHAPER™
VERSION	1.0
SHORT DESCRIPTION	Planning software for 2D/3D design and biofabrication process management
CONTENT	<p>License for SHAPER™ valid for 12 months (incl. patches updates, and upgrades), annually renewable</p> <p>Dongle-based license (USB) allowing remote and on-site design work</p> <p>Compatible with R-GEN instruments (R-GEN 100 and R-GEN 200)</p> <p>Additional software modules available for kits and tools</p>
PACKAGES	3, 5 or 10 Independent licences,
HARDWARE & OS REQUIREMENTS	
OPERATING SYSTEM	Windows 10 64-bit
SYSTEM RAM	8 Gb (minimum), 16 Gb (recommended)
HARD DISK STORAGE	512 Gb SSD (recommended)
CPU	<p>Minimum: Intel Core i7-4500U, 64-bit, 1.80 GHz (3.00 GHz Turbo), 4 MB cache, 2 cores, 4 threads</p> <p>Recommended: Intel Core i7-9850H, 64-bit, 2.60 GHz (4.60 GHz Turbo), 12 MB cache, 6 cores, 12 threads</p>
USB	1 free USB port (for the license dongle)
SUPPORTED 2D/3D FILE FORMAT	*.STL
ADDITIONAL HW	USB port available, mouse and keyboard

All the data/information are subject to change without notice